

IDA-1S One-Channel Infusion Device Analyzer

Technical Data



Ensure infusion pumps are tested accurately and quickly with real time results using the IDA-1S One-Channel Infusion Device Analyzer, the newest member of the Fluke Biomedical IDA family. The IDA-1S is a portable, battery-operated instrument that allows for speedy verification of infusion device performance.

The IDA-1S measures the flow rate and volume delivered, and the pressure generated in occlusions or blockages of the fluid line. The IDA-1S is based on sophisticated measurement technology trusted by biomedical professionals around the world. It is easy to set up and requires little or no training to use. The IDA-1S can be used to test a wide variety of infusion pumps and an auto-start feature simplifies syringe pump testing and other tests that have long startup times.

Key features:

- Integrated carrying handle and lightweight (2.7 lb) for easy portability
- Battery powered with up to 10 hours of continuous operation for on-the-go operation
- LCD touch screen for ease of use
- Average and instantaneous flow measurement
- Occlusion pressure measurements to 45 psi
- Maximize accuracy with Autostart mode enabling unit to begin testing only when fluid is detected
- Compatible with a wide variety of infusion pumps
- Based on technology that is proven and trusted worldwide
- On-board memory allows test results storage instantly
- Hydrograph graphical software to control unit, display results and print results via PC
- Global sales, service and support

Technical specifications

Flow rate measurement	
Technique	Flow is calculated by measuring volume over time
Range	0.5 ml/h to 1000 ml/h
Accuracy	1 % of reading \pm 1 LSD for flows of 16 ml/h to 200 ml/h for volumes over 20 ml; otherwise 2 % of reading \pm 1 LSD for volumes over 10 ml under laboratory conditions
Max test duration	10 hours on battery
Volume measurement	
Technique	Volume is measured directly by the measuring module in minimum sample sizes of 60 μ l
Range	0.06 ml to 999 ml
Accuracy	1 % of reading \pm 1 LSD for flow rates of 16 ml/h to 200 ml/h for volumes over 20 ml. Otherwise 2 % of reading \pm 1 LSD for volumes over 10 ml under laboratory conditions
Max test duration	10 hours on battery
Pressure measurement	
Technique (Occlusion test)	Direct measurement of pressure at the inlet port
Range	0 psi to 45 psi and equivalent in mmHg, Bar and kPa
Accuracy	1 % of Full Scale \pm 1 LSD under laboratory conditions
Max test duration	30 minutes
General specifications	
Storage of results	Test results stored for later viewing, printing or transfer to PC. Typical practical capacity: 100 tests
Power down	The results of tests in progress will be saved in the case of accidental power down
Computer control	The Product can be fully controlled from a PC using HydroGraph V3 software for IDA-1S
Battery power	4 x Panasonic HHR210AB NiMh 2000 mAh batteries
Charger	Operating voltage range: 100 V ac to 240 V ac
	Supply frequency: 50 Hz / 60 Hz
	Supply power: <20 VA
Size (HxWxD)	30 cm x 17 cm x 10 cm (12 in x 8 in x 4 in)
Weight	~1.2 kg (2.7 lb)
Temperature	Operating: 15 °C to 30 °C (59 °F to 86 °F)
	Storage: -20 °C to +40 °C (-4 °F to +104 °F) when drained of all liquid
Humidity	10 % to 90 % non-condensing
Altitude	0 meters to 2000 meters (6500 feet)
Safety	IEC 61010-1: Overvoltage category II, Pollution Degree 2
Electromagnetic environment	IEC 61326-1: Basic
Emissions classification	IEC CISPR 11: Group 1, Class A. Group 1 have intentionally generated and/or use conductively coupled radio-frequency energy which is necessary for the internal functioning of the equipment itself. Class A equipment is suitable for use in nondomestic locations and/or directly connected to a low-voltage power supply network
FCC	CFR47: Class A Part 15 subpart B
Electromagnetic compatibility	Applies to use in Korea only. Class A: Equipment (Industrial Broadcasting and Communication Equipment) ¹

This product meets requirements for industrial (Class A) electromagnetic wave equipment and the seller or user should take notice of it. This equipment is intended for use in business environments and is not to be used in homes.

Ordering information

Models/descriptions

4468525 IDA-1S One-Channel Infusion Device Analyzer

Standard accessories

4418071 Hydrograph Software and Users Manual
4497350 20 ml syringe
4480194 3-way plastic Luerlock
4478942 Drain tube (1 m)
4541948 Micro-90 bottle (225 ml)
1740487 USB data transfer cable
2461300 4-plug mains adapter kit for US, UK, EURO, AUS
4329971 Power supply cord without main adapters, univ, wall mount
4481150 NiMH replacement battery
2523995 Soft carrying case

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical.

Trusted for the measurements that matter.

Fluke Biomedical

6045 Cochran Road
Cleveland, OH 44139-3303 U.S.A.

For more information, contact us at:
(800) 850-4608 or Fax (440) 349-2307
Email: sales@flukebiomedical.com
Web access: www.flukebiomedical.com

©2014,2015 Fluke Biomedical.
Specifications subject to change without notice.
Printed in U.S.A. 9/2015 6002180c-en

**Modification of this document is not permitted
without written permission from Fluke Corporation.**