

**FLUKE**®

— **Biomedical**

# ProSim™ 8

Vital Signs Simulator

**Getting Started Manual**

PN 3984515

January 2011, Rev. 3, 3/16

© 2011-2016 Fluke Corporation. All rights reserved. Specifications are subject to change without notice.

All product names are trademarks of their respective companies.

## **Warranty and Product Support**

Fluke Biomedical warrants this instrument against defects in materials and workmanship for one year from the date of original purchase OR two years if at the end of your first year you send the instrument to a Fluke Biomedical service center for calibration. You will be charged our customary fee for such calibration. During the warranty period, we will repair or at our option replace, at no charge, a product that proves to be defective, provided you return the product, shipping prepaid, to Fluke Biomedical. This warranty covers the original purchaser only and is not transferable. The warranty does not apply if the product has been damaged by accident or misuse or has been serviced or modified by anyone other than an authorized Fluke Biomedical service facility. NO OTHER WARRANTIES, SUCH AS FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSED OR IMPLIED. FLUKE SHALL NOT BE LIABLE FOR ANY SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OR LOSSES, INCLUDING LOSS OF DATA, ARISING FROM ANY CAUSE OR THEORY.

This warranty covers only serialized products and their accessory items that bear a distinct serial number tag. Recalibration of instruments is not covered under the warranty.

This warranty gives you specific legal rights and you may also have other rights that vary in different jurisdictions. Since some jurisdictions do not allow the exclusion or limitation of an implied warranty or of incidental or consequential damages, this limitation of liability may not apply to you. If any provision of this warranty is held invalid or unenforceable by a court or other decision-maker of competent jurisdiction, such holding will not affect the validity or enforceability of any other provision.

# Notices

---

## All Rights Reserved

© Copyright 2016, Fluke Biomedical. No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language without the written permission of Fluke Biomedical.

---

## Copyright Release

Fluke Biomedical agrees to a limited copyright release that allows you to reproduce manuals and other printed materials for use in service training programs and other technical publications. If you would like other reproductions or distributions, submit a written request to Fluke Biomedical.

---

## Unpacking and Inspection

Follow standard receiving practices upon receipt of the instrument. Check the shipping carton for damage. If damage is found, stop unpacking the instrument. Notify the carrier and ask for an agent to be present while the instrument is unpacked. There are no special unpacking instructions, but be careful not to damage the instrument when unpacking it. Inspect the instrument for physical damage such as bent or broken parts, dents, or scratches.

---

## Technical Support

For application support or answers to technical questions, either email [techservices@flukebiomedical.com](mailto:techservices@flukebiomedical.com) or call 1-800- 850-4608 or 1-440-248-9300. In Europe, email [techsupport.emea@flukebiomedical.com](mailto:techsupport.emea@flukebiomedical.com) or call +31-40-2965314.

---

## Claims

Our routine method of shipment is via common carrier, FOB origin. Upon delivery, if physical damage is found, retain all packing materials in their original condition and contact the carrier immediately to file a claim. If the instrument is delivered in good physical condition but does not operate within specifications, or if there are any other problems not caused by shipping damage, please contact Fluke Biomedical or your local sales representative.

---

## Returns and Repairs

### Return Procedure

All items being returned (including all warranty-claim shipments) must be sent freight-prepaid to our factory location. When you return an instrument to Fluke Biomedical, we recommend using United Parcel Service, Federal Express, or Air Parcel Post. We also recommend that you insure your shipment for its actual replacement cost. Fluke Biomedical will not be responsible for lost shipments or instruments that are received in damaged condition due to improper packaging or handling.

Use the original carton and packaging material for shipment. If they are not available, we recommend the following guide for repackaging:

- Use a double-walled carton of sufficient strength for the weight being shipped.
- Use heavy paper or cardboard to protect all instrument surfaces. Use nonabrasive material around all projecting parts.
- Use at least four inches of tightly packed, industry-approved, shock-absorbent material around the instrument.

### Returns for partial refund/credit:

Every product returned for refund/credit must be accompanied by a Return Material Authorization (RMA) number, obtained from our Order Entry Group at 1-440-498-2560.

### Repair and calibration:

To find the nearest service center, go to [www.flukebiomedical.com/service](http://www.flukebiomedical.com/service) or

In the U.S.A. and Asia:

Cleveland Calibration Lab

Tel: 1-800-850-4608 x2564

Email: [globalcal@flukebiomedical.com](mailto:globalcal@flukebiomedical.com)

In Europe, Middle East, and Africa:

Eindhoven Calibration Lab

Tel: +31-40-2675300

Email: [ServiceDesk@fluke.com](mailto:ServiceDesk@fluke.com)

To ensure the accuracy of the Product is maintained at a high level, Fluke Biomedical recommends the product be calibrated at least once every 12 months. Calibration must be done by qualified personnel. Contact your local Fluke Biomedical representative for calibration.

---

## Certification

This instrument was thoroughly tested and inspected. It was found to meet Fluke Biomedical's manufacturing specifications when it was shipped from the factory. Calibration measurements are traceable to the National Institute of Standards and Technology (NIST). Devices for which there are no NIST calibration standards are measured against in-house performance standards using accepted test procedures.

---

## WARNING

Unauthorized user modifications or application beyond the published specifications may result in electrical shock hazards or improper operation. Fluke Biomedical will not be responsible for any injuries sustained due to unauthorized equipment modifications.

---

## Restrictions and Liabilities

Information in this document is subject to change and does not represent a commitment by Fluke Biomedical. Changes made to the information in this document will be incorporated in new editions of the publication. No responsibility is assumed by Fluke Biomedical for the use or reliability of software or equipment that is not supplied by Fluke Biomedical, or by its affiliated dealers.

---

## Manufacturing Location

The ProSim™ 8 Vital Signs Simulator is manufactured at Fluke Biomedical, 6920 Seaway Blvd., Everett, WA, U.S.A.



## ***Introduction***

The Fluke Biomedical ProSim™ 8 Vital Signs Simulator (hereafter the Product) is a full-featured, compact, portable simulator, used to measure the performance of patient monitors.

## ***Intended Use***

The Product is intended to be used to test and verify the basic operation of patient monitoring devices or systems used to monitor various physiological parameters of a patient, including ECG, Respiration, Invasive blood pressure, Non-invasive blood pressure, Temperature, and Cardiac output. Additionally, the Devices provide an optical signal to verify that the electronics within the pulse oximeter probe are functional.

The intended user is a trained biomedical equipment technician who performs periodic preventative maintenance checks on patient monitors in service. Users can be associated with hospitals, clinics, original equipment manufacturers and independent service companies that repair and service medical equipment. The end user is an individual, trained in medical instrumentation technology.

This Product is intended to be used in the laboratory environment, outside of the patient care area, and is not intended for use on patients, or to test devices while connected to patients. This Product is not intended to be used to calibrate medical equipment. It is intended for over the counter use.

## **Safety Information**

In this manual, a **Warning** identifies hazardous conditions and actions that could cause bodily harm or death. A **Caution** identifies conditions and actions that could damage the Product, the equipment under test, or cause permanent loss of data.

### **⚠⚠ Warnings**

**To prevent personal injury, use the Product only as specified, or the protection supplied by the Product can be compromised.**

**To prevent possible electrical shock, fire, or personal injury:**

- **Read all safety Information before you use the Product.**
- **Do not touch voltages >30 V ac rms, 42 V ac peak, or 60 V dc.**
- **Do not use the Product around explosive gas, vapor, or in damp or wet environments.**

- **Examine the case before you use the Product. Look for cracks or missing plastic. Carefully look at the insulation around the terminals.**
- **Do not use the Product if it is damaged.**
- **Disable the Product if it is damaged.**
- **Use this Product indoors only.**
- **Do not connect directly to mains.**
- **Use the correct terminals, function, and range for measurements.**
- **Do not use the Product if it operates incorrectly.**
- **Use only current probes, test leads, and adapters supplied with the Product.**
- **Remove all probes, test leads, and accessories that are not necessary for the measurement.**



- **Remove all probes, test leads, and accessories before the battery door is opened.**
- **The battery door must be closed and locked before you operate the Product.**
- **Replace the batteries when the low battery indicator shows to prevent incorrect measurements.**
- **Remove the batteries if the Product is not used for an extended period of time, or if stored in temperatures above 50 °C. If the batteries are not removed, battery leakage can damage the Product.**
- **Use only the mains power cord and connector approved for the voltage and plug configuration in your country and rated for the Product.**
- **Replace the mains power cord if the insulation is damaged or if the insulation shows signs of wear.**
- **Use only the external mains power supply included with the Product.**
- **Connect the battery charger to the mains power outlet before the Product.**
- **Do not put metal objects into connectors.**
- **Do not connect the Product to a patient or equipment connected to a patient. The Product is intended for equipment evaluation only and should never be used in diagnostics, treatment, or any other capacity where the Product would come in contact with a patient.**

- **Batteries contain hazardous chemicals that can cause burns or explode. If exposure to chemicals occurs, clean with water and get medical aid.**
- **Do not disassemble the battery.**
- **Do not disassemble or crush battery cells and battery packs.**
- **Do not put battery cells and battery packs near heat or fire. Do not put in sunlight.**
- **Do not short the battery terminals together.**
- **Do not keep cells or batteries in a container where the terminals can be shorted.**

- **Remove the input signals before you clean the Product.**
- **Use only specified replacement parts.**
- **Have an approved technician repair the Product.**

 **Warnings**












**For safe operation and maintenance of the Product:**

- **Repair the Product before use if the battery leaks.**
- **Use only Fluke Biomedical approved power adapters to charge the battery.**

## Symbols


Table 1 describes symbols used in association with the Product.

**Table 1. Symbols**

Symbol	Description	Symbol	Description
	WARNING - RISK OF DANGER. Consult user documentation.		WARNING. HAZARDOUS VOLTAGE. Risk of electric shock.
	Magnetic Field		Input jack for the DC output of the AC-DC supply connector.
	Conforms to European Union directives.		Certified by CSA Group to North American safety standards.
	Conforms to relevant Australian EMC standards.		Complies with 47 CFR Part 15 requirements of the U.S. Federal Communications Commission.
	Conforms to relevant South Korean EMC Standards.		Spent Lithium batteries should be disposed of by a qualified recycler or hazardous materials handler per local regulations. Contact your authorized Fluke Service Center for recycling information.
	This product complies with the WEEE Directive marking requirements. The affixed label indicates that you must not discard this electrical/electronic product in domestic household waste. Product Category: With reference to the equipment types in the WEEE Directive Annex I, this product is classed as category 9 "Monitoring and Control Instrumentation" product. Do not dispose of this product as unsorted municipal waste.		

## **Instrument Familiarization**

Table 2 is a list of Product top-panel controls and connections shown in Figure 1.



glh034.eps

**Figure 1. Top-Panel Controls and Connections**

**Table 2. Top-Panel Controls and Connections**

<b>Item</b>	<b>Name</b>	<b>Description</b>
①	ECG Posts	Connection posts for Device Under Test (DUT) ECG leads.
②	ECG Function	Accesses the ECG waveforms (adult, pediatric, and arrhythmias), and ECG test functions (performance waves, QRS detection, Tall T wave rejection, and R wave detection).
③	NIBP Button	Accesses the Non-Invasive Blood Pressure (NIBP) functions.
④	Special Functions	Accesses the temperature, respiration, cardiac output, fetal simulation, autosequences, and view memory functions.
⑤	SETUP Button	Accesses the setup controls.
⑥	IBP Button	Accesses the Invasive Blood Pressure (IBP) functions.
⑦	Navigation Buttons	Cursor control buttons for navigating menus and lists.
⑧	Enter Button	Sets the highlighted function.
⑨	Backlight Button	Turns the display backlight on and off.
⑩	Power Button	Turns the Product on and off.
⑪	SpO2 Button	Accesses the SpO2 functions.
⑫	Function Softkeys	Keys F1 through F5 are used to select from a number of selections that appear in the LCD display above each function softkey.

**Table 2. Top-Panel Controls and Connections (cont.)**

<b>Item</b>	<b>Name</b>	<b>Description</b>
⑬	LCD Display	Color display.
⑭	SpO2 Connector	Connector to the SpO2 accessory.
⑮	IBP Channel 2 Connector	Connector to an IBP input of the patient monitor.
⑯	Cardiac Output Connector	Connector to the Cardiac input of the patient monitor.
⑰	IBP Channel 1 Connector	Connector to the IBP input of the patient monitor.
⑱	Temperature Connector	Connector to the Temperature input of the patient monitor.

Table 3 is a list of Product controls and connections shown in Figure 2.




Figure 2. Back, Front, and Side Panel Connections

glh035.eps

**Table 3. Back, Front, and Side Panel Connections**

<b>Item</b>	<b>Name</b>	<b>Description</b>
①	AC/DC Supply Connector	Input jack for the DC output of the AC/DC supply connector.
②	Battery Charge LED	Battery charges when LED shows red. Green shows battery charge is complete.
③	Battery Latch	Locks the battery pack into the Product. Push down to remove the battery pack.
④	Mini B USB Device Port	Used to connect to a PC for remote control or download test results data to a PC.
⑤	USB A Controller Port	For external keyboard, barcode reader, or printer.
⑥	ECG BNC Connector	High-level output of ECG signal.
⑦	Air Port Connector	Pressure port for NIBP cuff and monitor.
⑧	Magnetic Holder for SpO2 Finger Module	Holds the SpO2 Optical Emitter and Detector finger module in two orientations.



## **How to Turn the Product On**

After you unpack and inspect the Product, fully charge the battery before the first use. Afterwards, charge the battery when the Product shows the low battery message. See the Users Manual for detailed instructions.

Push **Ⓞ** on the front panel to turn the Product on. The startup screen shows on the display. When the self test is complete and no errors are sensed, the Home screen shows in the display. The Product is ready to use. See the Users Manual for detailed instructions.

## **General Specifications**

### **Temperature**

Operating ..... 10 °C to 40 °C (50 °F to 104 °F)  
Storage..... -20 °C to +60 °C (-4 °F to +140 °F)

**Humidity** ..... 10 % to 90 % non-condensing

**Altitude** ..... 3000 meters (9843 ft)

**Size (L x W x H)**..... 30.22 cm x 14.48 cm x 8.64 cm (11.9 in x 5.7 in x 3.4 in)

**Display**..... LCD Color Display

### **Communication**

USB Device Upstream Port..... Mini-B connector for control by a computer

USB Host Controller Port ..... Type A, 5 V output, 0.5 A max load. Connector for keyboard, barcode reader, and printer

Wireless..... IEEE 802.15.4 for control by a computer

#### USB Device Virtual COM Port Settings

Baud Rate ..... 115 200 bps

Data bits ..... 8 data bits

Stop Bits..... 1 stop bit

Flow Control ..... Hardware (RTS/CTS)

<b>Power</b> .....	Lithium-Ion rechargeable, 7.2 V, 31 Wh battery, 4300 mAh
<b>Battery Charger</b> .....	100 V to 240 V, 50/60 Hz input, 15 V/2.0 A output. For best performance, the battery charger should be connected to a properly grounded ac receptacle
<b>Battery Life</b> .....	9 hours (minimum), 100 NIBP cycles typical
<b>Weight</b> .....	1.81 kg (4 lb)

**Wireless Radio**

Frequency Range .....	2412 MHz to 2462 MHz
Output Power .....	<1 mW

**Safety** ..... IEC 61010-1: Overvoltage Category II, Pollution Degree 2

**Electromagnetic Compatibility (EMC)**

International ..... IEC 61326-1: Portable Electromagnetic Environment CISPR 11: Group 1, Class A

*Group 1: Equipment has intentionally generated and/or uses conductively-coupled radio frequency energy that is necessary for the internal function of the equipment itself.*

*Class A: Equipment is suitable for use in all establishments other than domestic and those directly connected to a low-voltage power supply network that supplies buildings used for domestic purposes. There may be potential difficulties in ensuring electromagnetic compatibility in other environments due to conducted and radiated disturbances.*

*Emissions that exceed the levels required by CISPR 11 can occur when the equipment is connected to a test object. The equipment may not meet the immunity requirements of this standard when test leads and/or test probes are connected.*

Korea (KCC) ..... Class A Equipment (Industrial Broadcasting & Communication Equipment)

*Class A: Equipment meets requirements for industrial electromagnetic wave equipment and the seller or user should take notice of it. This equipment is intended for use in business environments and not to be used in homes.*

USA (FCC) .....47 CFR 15 Intentional Radiators: This device complies with part 15 of the FCC Rules.

*Operation is subject to the following two conditions:*

*(1) This device may not cause harmful interference.*

*(2) This device must accept any interference received, including interference that may cause undesired operation. (15.19).*

*Changes or modifications not expressly approved by Fluke could void the user's authority to operate the equipment. (15.21)*

## Detailed Specifications

### Normal-Sinus-Rhythm Waveform

**ECG Reference** .....The ECG amplitudes specified are for Lead II (calibration), from the baseline to the peak of the R wave. All other leads are proportional.

**Normal Sinus Rhythm** .....12-lead configuration with independent outputs referenced to right leg (RL). Output to 10 Universal ECG Jacks, color-coded to AHA and IEC Standards.

**High-Level Output** .....0.5 V/mV  $\pm$ 5 % of the ECG amplitude setting available on a BNC connector.

**Amplitude** .....0.05 mV to 0.5 mV (0.05 mV steps); 0.5 mV to 5.0 mV (0.25 mV steps). Other leads are proportional to Lead II (reference lead) in percentage per:

Lead I	70
Lead II	100
Lead III	30
Lead V1	24
Lead V2	48
Lead V3	100
Lead V4	120
Lead V5	112
Lead V6	80

<b>Amplitude Accuracy</b> .....	±(2 % of setting + 0.05 mV)
<b>ECG Rate</b> .....	10 BPM to 360 BPM in 1 BPM steps
<b>Rate Accuracy</b> .....	±1 % of setting
<b>ECG Waveform Selection</b> .....	Adult (80 ms) or pediatric (40 ms) QRS duration
<b>ST-Segment Elevation</b> .....	Adult mode only. -0.8 mV to +0.8 mV (0.1 mV steps) Additional steps: +0.05 mV and -0.05 mV
<b>Power-On Default</b> .....	60 BPM, 1.0 mV, adult QRS and ST-segment elevation of 0 mV

***Pacemaker Waveform***

**Pacer-Pulse**

<b>Amplitude</b> .....	0 mV (off), ±2 mV, ±4 mV, ±6 mV, ±8 mV, ±10 mV, ±12 mV, ±14 mV, ±16 mV, ±18 mV, ±20 mV, ±50 mV, ±100 mV, ±200 mV, ±500 mV, and ±700 mV for lead II (reference lead)
------------------------	--

**Accuracy**

Reference lead II .....	±(5 % setting + 0.2 mV)
All other leads .....	±(10 % setting + 0.4 mV)

<b>Pacer-Pulse Width</b> .....	0.1 ms, 0.2 ms, 0.5 ms, 1.0 ms, and 2.0 ms ±5 %
--------------------------------	---

<b>Paced Arrhythmias</b> .....	Atrial 80 BPM Asynchronous 75 BPM Demand with frequent sinus beats Demand with occasional sinus beats Atrio-Ventricular sequential Noncapture (one time) Nonfunction
--------------------------------	--

<b>Power-On Default</b> .....	Amplitude 10 mV, width 1.0 ms, atrial waveform
-------------------------------	--

**Arrhythmia**

<b>Baseline NSR</b> .....	80 BPM
<b>PVC Focus</b> .....	Left focus, standard timing (except where specified)
<b>Supraventricular Arrhythmia</b> .....	Atrial fibrillation (coarse or fine); atrial flutter; sinus arrhythmia; missed beat (one time); atrial tachycardia; paroxysmal atrial tachycardia; nodal rhythm; and supraventricular tachycardia.
<b>Premature Arrhythmia</b> .....	Premature atrial contraction (PAC); premature nodal contraction (PNC); PVC1 left ventricular; PVC1 left ventricular, early; PVC1 left ventricular, R on T; PVC2 right ventricular; PVC2 right ventricular, early; PVC2 right ventricular, R on T; and multifocal PVCs
<b>Ventricular Arrhythmia</b> .....	PVCs 6 per minute, 12 per minute, or 24 per minute; frequent multifocal PVCs; bigeminy; trigeminy; multiple PVCs (one-time run of 2 PVCs, 5 PVCs, or 11 PVCs); mono-ventricular tachycardia (120 BPM to 300 BPM in 5 BPM steps); poly-ventricular tachycardia (5 types); ventricular fibrillation (coarse or fine); and Asystole
<b>Conduction Defect</b> .....	First-, second-, or third-degree heart block; and right- or left-bundle-branch block
<b>Advanced Cardiac Life Support</b>	
Shockable Pulseless Arrest Rhythms .....	Ventricular fibrillation (coarse), ventricular fibrillation (fine), unstable polymorphic ventricular tachycardia
Non Shockable Pulseless Arrest Rhythms.....	Asystole
Symptomatic Bradycardia .....	Sinus Bradycardia (<60 BPM) 2 <sup>nd</sup> Degree AV Block, Mobitz Type I 2 <sup>nd</sup> Degree AV Block, Mobitz Type II Complete/3 <sup>rd</sup> Degree AV Block Right Bundle Branch Block Left Bundle Branch Block

Symptomatic Tachycardia

Regular Narrow-complex Tachycardias (QRS <0.12 seconds)

Sinus Tachycardia ..... >150 BPM

Supraventricular Tachycardia ..... SVT

Regular Wide-complex Tachycardias (QRS ≥0.12 seconds)

Sinus Tachycardia ..... >150 BPM

Supraventricular Tachycardia ..... SVT with aberrancy

Irregular Tachycardia ..... Atrial Fibrillation (Coarse and fine), Atrial Flutter, unstable monomorphic ventricular tachycardia (120 BPM to 300 BPM), Torsade De Pointes/Polymorphic ventricular tachycardia (long QT interval)

**ECG-Performance-Testing**

**Amplitude (peak-to-peak)** ..... 0.05 mV to 0.5 mV (0.05 mV steps) 0.5 mV to 5.0 mV (0.25 mV steps) Other leads are proportional to Lead II (reference lead) in percentage per:

Lead I ..... 70

Lead II ..... 100

Lead III ..... 30

Lead V1 through V6 ..... 100

**Pulse Wave** ..... 30 BPM, 60 BPM, with 60 ms pulse width

**Square Wave** ..... 0.125 Hz, 2.0 Hz, 2.5 Hz

**Triangle Wave** ..... 0.125 Hz, 2.0 Hz, 2.5 Hz

**Sine Wave** ..... 0.05 Hz, 0.5 Hz, 1 Hz, 2 Hz, 5 Hz, 10 Hz, 25 Hz, 30 Hz, 40 Hz, 50 Hz, 60 Hz, 100 Hz, and 150 Hz

**R-wave Detection**

Waveform .....Triangular pulse  
Rate.....30 BPM, 60 BPM, 80 BPM, 120 BPM, 200 BPM, and 250 BPM  
Width .....8 ms to 20 ms in 2 ms steps, and 20 ms to 200 ms in 10 ms steps  
Width Accuracy ..... $\pm(1\% \text{ of setting} + 1 \text{ ms})$

**QRS Detection**

Widths .....8 ms to 20 ms in 2 ms steps and 20 ms to 200 ms in 10 ms steps  
Width Accuracy ..... $\pm(1\% \text{ of setting} + 1 \text{ ms})$   
Rate.....30 BPM, 60 BPM, 80 BPM, 120 BPM, 200 BPM, and 250 BPM  
R-Wave up slope.....0.875 amplitude, 0.4375 x width  
R-Wave down slope .....Full amplitude, 0.5 x width  
S-Wave up slope.....0.125 amplitude, 0.0625 x width

**Tall T-Wave Rejection**

Waveform  
    QT Interval .....350 ms  
    T-Wave width .....180 ms  
    T-Wave Shape ..... $\frac{1}{2}$  sinewave  
Amplitude .....0 % to 150 % reference lead amplitude in 10 % steps  
Rate.....80 BPM

**Rate Accuracy** ..... $\pm 1\% \text{ of setting}$

**Amplitude Accuracy**..... $\pm(2\% \text{ of setting} + 0.05 \text{ mV})$

### **ECG Artifact**

Type.....	50 Hz, 60 Hz, muscular, baseline wander, respiration
Size.....	25 %, 50 %, 100 % of the normal sinus R-Wave for each lead
Lead Select.....	All, RA, LL, LA, V1, V2, V3, V4, V5, V6

### **Fetal / Maternal-ECG**

Fetal Heart Rate (Fixed).....	60 BPM to 240 BPM in 1 BPM steps
Fetal Heart Rate (IUP).....	140 BPM at beginning, then varies with pressure
Intrauterine-Pressure Waveforms.....	Early deceleration, late deceleration, and acceleration
Wave Duration.....	90 seconds, bell-shaped pressure curve, from 0 mmHg to 90 mmHg and returning to 0 mmHg
IUP Period.....	2 minutes, 3 minutes, or 5 minutes; and manual
Default Settings.....	FHR 140 BPM, early deceleration wave, manual

### **Invasive Blood Pressure**

Channels.....	2, each independently settable with identical parameters and are individually electrically isolated from all other signals
Input/output Impedance.....	300 $\Omega$ $\pm$ 10 %
Exciter Input Range.....	2.0 V to 16.0 V peak
Exciter-Input Frequency Range.....	DC to 5000 Hz
Transducer Sensitivity.....	5 $\mu$ V/V/mmHg (default) or 40 $\mu$ V/V/mmHg
Pressure Accuracy.....	$\pm$ (1 % of setting + 1 mmHg) Accuracy guaranteed for DC excitation only
Static Pressure.....	-10 mmHg to +300 mmHg in 1 mmHg steps
Pressure Units.....	mmHg or Kpa



**Dynamic Waveforms**

Types (default pressures) ..... Arterial (120/80)  
 Radial artery (120/80)  
 Left ventricle (120/00)  
 Right ventricle (25/00)  
 Pulmonary artery (25/10)  
 Pulmonary-artery wedge (10/2)  
 Right atrium (central venous or CVP) (15/10)

Pressure Variability ..... Systolic and diastolic pressures are independently variable in 1 mmHg steps.

**Swan-Ganz Sequence** ..... Right atrium, right ventricle (RV), pulmonary artery (PA), pulmonary artery wedge (PAW)

**Cardiac Catheterization**

Chambers ..... Aortic, Pulmonary valve, and Mitral valve

**Respiration Artifact**

Arterial, radial artery, and left ventricle ..... 5 % to 10 % multiplication

Other ..... 5 mmHg or 10 mmHg

**BP Output** ..... Circular DIN 5-Pin

**Power-On Default** ..... 0 mmHg

**Respiration**

**Rate** ..... 0 BrPM (OFF), 10 BrPM to 150 BrPM in 1 BrPM steps

**Waves** ..... Normal or ventilated

**Ratio (inspiration:expiration)**

Normal ..... 1:1, 1:2, 1:3, 1:4, 1:5

Ventilated ..... 1:1

**Impedance Variations ( $\Delta \Omega$ )** ..... 0.00  $\Omega$  to 1.00  $\Omega$  in 0.05  $\Omega$  steps and 1.00  $\Omega$  to 5.00  $\Omega$  in 0.25  $\Omega$  steps

**Delta Accuracy** .....  $\pm$ (5 % of setting + 0.1  $\Omega$ )

<b>Baseline</b> .....	500 $\Omega$ , 1000 $\Omega$ (default), 1500 $\Omega$ , 2000 $\Omega$ , Leads I, II, III
<b>Baseline Accuracy</b> .....	$\pm 5$ %
<b>Respiration Lead</b> .....	LA or LL (default)
<b>Apnea Selection</b> .....	12 seconds, 22 seconds, or 32 seconds (one-time events), or continuous (Apnea ON = respiration OFF)
<b>Power-On Default</b> .....	20 BrPM, delta 1.0 $\Omega$

### **Temperature**

<b>Temperature</b> .....	30.0 $^{\circ}\text{C}$ to 42.0 $^{\circ}\text{C}$ in 0.5 $^{\circ}\text{C}$ steps
<b>Accuracy</b> .....	$\pm 0.4$ $^{\circ}\text{C}$
<b>Compatibility</b> .....	Yellow Springs, Inc. (YSI) Series 400 and 700
<b>Output</b> .....	Circular DIN 4-pin

### **Cardiac Output**

<b>Catheter Type</b> .....	Baxter Edwards, 93a-131-7f
<b>Calibration Coefficient</b> .....	0.542 (0 $^{\circ}\text{C}$ injectate), 0.595 (24 $^{\circ}\text{C}$ injectate)
<b>Blood Temperature</b> .....	36 $^{\circ}\text{C}$ (98.6 $^{\circ}\text{F}$ ) to 38 $^{\circ}\text{C}$ (100.4 $^{\circ}\text{F}$ ) $\pm 0.2$ $^{\circ}\text{C}$ in 1 $^{\circ}\text{C}$ steps
<b>Injectate Volume</b> .....	10 cc
<b>Injectate Temperature</b> .....	0 $^{\circ}\text{C}$ or 24 $^{\circ}\text{C}$
<b>Cardiac Output</b> .....	2.5 liters per minute, 5 liters per minute, 10 liters per minute $\pm 7.5$ %
<b>Faulty-Injectate Curve</b> .....	Waveform for simulation available
<b>Left-to-Right-Shunt Curve</b> .....	Waveform for simulation available
<b>Calibrated Pulse</b> .....	1.5 $^{\circ}\text{C}$ for 1 second
<b>Connector</b> .....	Circular DIN 7 pin
<b>Power-On Default</b> .....	5 liters per minute, 0 $^{\circ}\text{C}$ injectate, 37 $^{\circ}\text{C}$ blood temperature

**Non-Invasive Blood Pressure**

**Pressure Units** .....mmHg or kPa

**Manometer (Pressure Meter)**

Range.....0 mmHg to 400 mmHg

Resolution .....0.1 mmHg

Accuracy..... $\pm(0.5\% \text{ reading} + 0.5 \text{ mmHg})$

**Pressure Source**

Target Pressure Range ..... 15 mmHg to 400 mmHg

Resolution ..... 1 mmHg

**NIBP Simulations**

Pulse .....2 mmHg max into 500 ml NIBP system

Volume of air moved ..... 1.25 ml max

Simulations.....Systolic/diastolic (MAP)

Adult .....60/30 (40), 80/50 (60), 100/65 (77); 120/80 (93); 150/100 (117); and 200/150 (167) and 255/195 (215)

Neonatal.....35/15 (22); 60/30 (40); 80/50 (60);100/65 (77);120/80 (93) and 150/100 (117)

Pressure variability.....Systolic and diastolic pressures are variable by 1 mmHg

Repeatability ..... Within  $\pm 2$  mmHg (at maximum pulse size independent of device under test)

**Synchronization**

Normal Sinus heart rates .....30 BPM to 240 BPM

Maximum rate at 1 ml.....240 BPM achievable with pulses up to 1 ml

Maximum rate at 1.25 ml.....180 BPM

Arrhythmias .....Premature atrial contraction (PAC), Premature ventricular contraction (PVC), atrial fibrillation, and missed beat.

**Leak Test**

Target Pressure .....	20 mmHg to 400 mmHg
Elapsed time .....	0:30 minutes:seconds to 5:00 minutes:seconds in 30 second steps
Range .....	0 mmHg/minute to 200 mmHg/minute
Internal Leak rate .....	<2 mmHg/min into 500 ml rigid volume
<b>Pressure Relief Test Range</b> .....	100 mmHg to 400 mmHg

***Oximeter SpO2 Optical Emitter and Detector (optional)***

**% O<sub>2</sub>**

Range .....	30 % to 100 %
Resolution .....	1 %

**Accuracy**

With oximeter manufacturer's R-curve

Saturation within UUT specific range.....	±(1 count + specified accuracy of the UUT)
Saturation outside UUT specific range .....	monotonic with unspecified accuracy

With Fluke Biomedical R-curves

91 % to 100 % .....	±(3 counts + specified accuracy of the UUT)
81 % to 90 % .....	±(5 counts + specified accuracy of the UUT)
71 % to 80 % .....	±(7 counts + specified accuracy of the UUT)
Below 71 % .....	monotonic with unspecified accuracy

**Heart Rate**

Range .....	30 BPM to 300 BPM in 1 BPM steps. Oximeter SpO2 Optical Emitter and Detector is synchronized with ECG rate delayed by 150 ms.
Accuracy .....	±1 % of setting

**Transmission (Ratio of detector current to LED current, expressed in parts per million (ppm))**

Range.....	0 ppm to 300.00 ppm
Resolution .....	0.01 ppm
Accuracy.....	+50 %/-30 % for compatible monitors, unspecified for others. Selected by finger size and color: Dark, thick finger, medium finger, light, thin finger, neonatal foot.

**Pulse Amplitude**

Range.....	0 % to 20.00 %
Resolution .....	0.01 %

**Artifact**

Respiration

Range.....	0 % to 5 % of transmission
Resolution .....	1 %
Rate.....	all ProSim respiration simulation settings

Ambient Light

Range.....	0X to 5X transmitted light
Resolution .....	1X
Frequency .....	DC, 50 Hz, 60 Hz, and 1 kHz to 10 kHz in 1 kHz steps

**Masimo Rainbow Technology**.....Masimo Rainbow technology with an optional adapter supplied by Masimo that allows the ProSim two wavelength to test the Rainbow multiple wavelength system

**Compatible Manufacturer Products**

With manufacturer R-curve .....	Nellcor, Masimo, Nonin, and Nihon Kohden
With Fluke Biomedical R-curve .....	Mindray, GE-Ohmeda, Philips/HP, and BCI

***Pre-Defined Simulations***

Normal  
Hypertensive  
Hypotensive  
Tachycardic  
Bradycardic  
Ventricular Fibrillation  
Asystole

***Autosequences (default)***

Monitor testing sequence  
Medical training sequence  
Oximeter testing sequence  
Cardiac failure sequence  
Arrhythmia sequence  
Exercise sequence  
Respiration sequence  
Performance Wave Test  
IBP testing sequence  
Temperature sequence