EasyLung

Instructions for use
Gebrauchsanweisung
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso

IMT Analytics


IMT Analytics AG · Gewerbestrasse 8 · 9470 Buchs · Switzerland
T +41 81 750 67 10 · www.imtanalytics.com

700.282.000_04, 2023-03

A

B

English - EN

Intended use

The EasyLung is intended to be used exclusively for regular functional testing of ventilators or anaesthetic machines.

The EasyLung is intended for exclusive use in manufacturing facilities, hospitals or other healthcare facilities.

The Adult silicone bag and the resistance head can be autoclaved at 134°C and 2 bar for a duration of 3.5 min as well as a drying time of at least 2 min. To avoid colour changes to the side panel, we recommend removing the EasyLung side panel beforehand. All parts of the test lung can be cleaned and disinfected with products suitable for plastics; we recommend using Gigasept FF solution. Before reusing, the EasyLung must be tested for function. To maintain the full functionality of the EasyLung, a tube system with a taper of 15 mm or 22 mm diameter should be used.

Contraindications

The test lung is intended exclusively for use with ventilators and anaesthetic machines. Direct use on patients is not permitted. Therefore, no contraindications are defined.

Specifications

Volume of the silicone bag	max. 1000 mL	Temperature range for operation: for transport and storage:	5°C–40°C –20°C–50°C
Compliance	25 mL/mbar, ± 15%	Humidity range for operation: for transport and storage:	10% RH–95% RH 10% RH–95% RH
Resistance	20 mbar/L/s; ± 15%	Atmospheric pressure range for operation: for transport and storage:	600 hPa–1100 hPa 500 hPa–1100 hPa
Reprocessing cycles	max. 100		

Free from latex, DEHP and BPA

A

Scope of delivery/Replacement parts

a. EasyLung side panel	700.064.000
b. Retainer clip	300.318.000
c. Adult silicone bag 1 L	301.185.000
d. Resistance head	300.759.000
e. Tension nut for the bag	300.761.000
Instructions for use	700.282.000
EasyLung complete	300.756.100

NOTE: Only the parts in the scope of delivery may be used. If the EasyLung needs to be transported, the original packaging must be used.

B

Disassembly/assembly

- Remove the retainer clip.
- Push the resistance head downwards.
- Pull the silicone bag out backwards and remove the tension nut.
- Pull the resistance head out sideways and remove the whole unit.
- The silicone bag can be detached from the resistance head.
- The lung is assembled in the reverse order. Finally, push the resistance head forwards and insert the retainer clip.
- To attach the lung to the machine, the EasyLung is positioned vertically on the mounting device.
- To attach the tubing system, the taper (22 mm or 15 mm) is placed on the resistance head.

CE conformity

This test lung satisfies the EU Medical Device Regulation (2017/745/EU).

Basic UDI-DI

764021464T01FP

Explanation of symbols

	Manufacturer's address
	LOT number
	Catalogue number

	European Representative: IMT Analytics GmbH, Fürst-Franz-Josef-Strasse 33, 9493 Mauren, Liechtenstein
	UK Responsible Person: Mr. Friberg 104B Elers road, London, W13 9QE England United Kingdom
	MDR Importer: IMT Analytics GmbH, Fürst-Franz-Josef-Strasse 33, 9493 Mauren, Liechtenstein
	Unique device identification
	CE marking
	Mark for United Kingdom Conformity Assessment
	Consult instructions for use
	Medical device
	Keep away from sunlight
	Comply with temperature range for transport and storage
	Comply with atmospheric humidity range for transport and storage
	Comply with air pressure range for transport and storage
	WEEE marking Waste Electrical and Electronic Equipment

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