

Technical data

INCU II™ Incubator/Radiant Warmer Analyzer

Protect a baby's safety by regularly performing all necessary preventive maintenance tests on infant incubators and radiant warmers. Life matters. Testing matters.

The Fluke Biomedical INCU II is an all-in-one Incubator/Radiant Warmer Analyzer that simplifies testing and ensures proper performance and safety of infant incubators, transport incubators and radiant warmers by complying to global standards (IEC 60601-2-19, IEC 60601-2-20, IEC 60601-2-21).

Portable and intuitive-to-use, the INCU II simultaneously tests temperature, airflow, sound, and humidity. Test results and parameters can be monitored in real-time on a large easy-to-read LCD screen, or be transferred to your computer wirelessly*.

Pass/Fail indicators are displayed at the end of tests for fast and easy troubleshooting, and a general testing option allows you to complete a quick test in just 15 minutes.

The INCU II can stay in the incubator chamber up to 24 hours using its rechargeable battery, or for 48 hours on AC power, without compromising the integrity of the test environment.



Key features

- Enables simultaneous measurement of humidity, airflow, sound, and temperature from six different points
- Displays pass/fail indicators and real-time test results on a large LCD screen
- Wireless* functionality supports fast and convenient downloading of test results and data
- Simplifies test set-up with color-coded temperature probes that match color of input ports
- Has both Standards-based testing and a General Testing option—which allows you to complete a test in as little as 15 minutes
- Portable and compact design, weighs less than three pounds
- Creates personalized test sequences for automatic test completion
- Able to choose from 10 different language interfaces
- Complies with all global standards, including IEC 60601-2-19, IEC 60601-2-20 and IEC 60601-2-21

Specifications

Physical specifications	
Size (HxWxL without sensors)	23 cm x 21 cm x 6 cm (9.0 in x 8.5 in x 2.5 in)
Total weight	12.2 cm x 9.7 cm x 4.8 cm (4.8 in x 3.8 in x 1.9 in)
Weight with sensors	1.4 kg (3 lb)
Pucks (5)	2.5 kg (5.5 lb)
Carrying case	1.1 kg (2.5 lb)
Power	
Power Adapter—Universal voltage	
Input	100 V to 240 V with adaptors 50 Hz/60 Hz
Output	15V dc, 1.3 A maximum
Rechargeable lithium-ion battery, internal	7.4 V, 7.8 Ah, 58 Wh 24 hour battery life with 30 second sample rate
Measurements and tests specifications	
Air convection peripheral temperature sensors for incubator (T1-T5)	
Number of sensors	5
Range	0 °C to 50 °C
Accuracy	±0.05 °C
Display resolution	0.01 °C
Relative humidity	
Range	0 % to 100 %
Accuracy	+/- 3 % RH (0 % to 100 % non-condensing)
Display resolution	0.1 % RH
Airflow	
Range	0.2 m/sec to 2.0 m/sec at 35°C, 50 % RH
Accuracy	+0.1 m/sec
Display resolution	0.01 m/sec
Sound pressure	
Range	30 dB(A) to 100 dB(A)
Accuracy	± 5 dB(A)
Display resolution	0.1 dB(A)
Standard	IEC-61672-1 Class 2 from 31.5Hz to 8kHz
Surface temperature	
Range	-5 °C to 60 °C
Accuracy	± 0.5 °C
Display resolution	0.05 °C

Specifications continued

Environmental specifications	
Operating temperature	10 °C to 40 °C
Storage temperature	-20 °C to 60 °C
Humidity	10 % to 90 % non-condensing
Altitude	2000 m
Ingress protection rating	IP-20
Safety	
IEC 61010-1	Overvoltage category none, pollution degree 2
IEC 61326-1	Basic
Electromagnetic compatibility (EMC)	
Emissions classification IEC CISPR11: Group 1, Class A.	Group 1 has intentionally generated and/or use conductively coupled radio-frequency energy, which is necessary for the internal functioning of the equipment itself. Class A equipment is suitable for use in nondomestic locations and/or directly connected to a low voltage power supply network.
USA (FCC) Intentional radiators	This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. (15.19)
Korea (KCC) Class A Equipment (Industrial Broadcasting & Communication Equipment)	Class A: Equipment meets requirements for industrial electromagnetic wave equipment and the seller or user should take notice of it. This equipment is intended for use in business environments and not to be used in homes.
Wireless module listing	
FCC (United States) compliant (Class A)	FCC ID: X3ZBTMOD3
IC (Industry Canada) compliant	IC: 8828A-MOD3
CE (European) certified	CE0051
802.15.1 qualified	QD ID: B019224
Wireless radio	
Frequency range	2412 MHz to 2483 MHz
Output power	10 mW
Interface	
Push navigation buttons	Power on/off, test, select, back, and arrow keys
User preferences	Adjust backlight, display brightness, and set clock
Recall and run templates on tester	General and personalized test templates allow for on-board test automation
Templates	Time duration, frequency of data capture and tests
Select user preferences	Units of measure, view test results of current and past tests on the Analyzer
View battery life	Indicator bar shows life remaining
Interface languages	English, French, German, Italian, Spanish, Brazilian Portuguese, Simplified Chinese, Japanese, Russian, Polish

Ordering information

Item #	Model	Description
4722620	INCU II-BT	Incubator Analyzer, Wireless
4722635	INCU II-NO BT	Incubator Analyzer, Non-Wireless

Standard accessories

4720996	INCU II-Cable	K-Type Thermocouple
4721000	INCU II S Sensor	Sound Sensor
4721017	INCU II A F Sensor	Air Flow Sensor
4721021	INCU II H Sensor	Humidity Sensor
4721039	INCU II R Sensor	Red Temperature Sensor
4721042	INCU II G Sensor	Green Temperature Sensor
4721056	INCU II Y Sensor	Yellow Temperature Sensor
4721063	INCU II W Sensor	White Temperature Sensor
4721074	INCU II B Sensor	Blue Temperature Sensor
4721109	INCU II Spacers	Tripod Spacers
4721111	INCU II R Puck	Red Temperature Puck
4721127	INCU II G Puck	Green Temperature Puck
4721130	INCU II Y Puck	Yellow Temperature Puck
4721148	INCU II W Puck	White Temperature Puck
4721153	INCU II B Puck	Blue Temperature Puck
4721166	INCU II USB	Type A to Micro B, USB Cable
4721194	INCU II Adapter	AC Adapter
4715713	INCU II Pad	Placement Pad Microfiber
4715749	INCU II Carry Case	Soft Carry Case
4724692	INCU II Puck Case	Soft Carry Case Puck

Optional accessories

4721175	INCU II ST Heater	Skin Temperature Heater Assembly
---------	-------------------	----------------------------------

Fluke Biomedical.

Trusted for the measurements that matter.

Fluke Biomedical
28775 Aurora Road
Cleveland, OH 44139 U.S.A.

For more information, contact us at:
(800) 850-4608 or Fax (440) 349-2307
Email: sales@flukebiomedical.com
Web access: www.flukebiomedical.com

©2015, 2017, 2019 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A.
12/2019 6006214e-en

Modification of this document is not permitted without written permission from Fluke Corporation.