

Mammography coverage table

Mammography machines supported by RaySafe, based on instrument models and parameters measured.



Parameters:

○ Dose, HVL & time

● kVp, dose, HVL & time

		Solo MAM Dose	Solo MAM	Xi MAM	Xi MAM M-Pro	Xi MAM M-Pro Plus	Xi Scanning MAM	X2 MAM Dose	X2 MAM	
Fujifilm	Amulet f/s/Innovality, Aspire HD-s/Cristalle/HD Plus	W/Rh	○	●	○	●	●	○	●	
		W/Al					●	○	●	
	Amulet, Aspire HD	Mo/Mo	○	●	●	●	●	○	●	
		Mo/Rh	○	○	○	●	●	○	●	
		W/Rh	○	●	○	●	●	○	●	
GE	Performa	Mo/Mo	○	●	●	●	●	○	●	
		Mo/Rh	○	○	○	●	●	○	●	
	Senographe Pristina	Mo/Mo	○	○	○	○	○	○	○	
		Rh/Ag							○	●
		Mo/Cu							○	○
	Senographe DS, Senographe Essential	Mo/Mo	○	●	●	●	●	○	○	
		Mo/Rh	○	○	○	●	●	○	○	
		Rh/Rh	○	○	○	○	○	○	○	
	Senographe DMR	Mo/Mo	○	●	●	●	●	○	○	
		Mo/Rh	○	○	○	●	●	○	○	
Mo/Al		○	○	○	○	○	○	○		
Rh/Rh		○	○	○	○	○	○	○		
	Rh/Al	○	○	○	○	○	○	○		
SenoBright	Mo/Mo	○	●	●	●	●	○	○		
	Mo/Rh	○	○	○	●	●	○	○		
	Rh/Rh	○	○	○	○	○	○	○		
	Mo/Cu						○	○		
	Rh/Cu						○	○		
Senographe Crystal	W/Rh	○	●	○	●	●	○	●		
Hologic	Selenia Dimensions	W/Rh	○	●	○	●	●	○	●	
		W/Ag	○	○		○	○	○	●	
		W/Al					●	○	●	
		Combo					○	○	●	
	Selenia (W)	W/Rh	○	●	○	●	●	○	●	
	W/Ag	○	○		○	○	○	●		
M-IV, Selenia (Mo)	Mo/Mo	○	●	●	●	●	○	●		
	Mo/Rh	○	○	○	●	●	○	●		
MultiCare	Mo/Mo	○	●	●	●	●	○	●		
IMS	Giotto 3D/3DL, Giotto Tomo	W/Rh	○	●	○	●	●	○	●	
		W/Ag	○	○		○	○	○	●	
	Giotto Class	W/Al						○	●	
		W/Ag						○	●	
		W/Cu						○	○	
Giotto MD/SD/SDL	W/Rh	○	●	○	●	●	○	●		
Giotto Image	Mo/Mo	○	●	●	●	●	○	●		
	Mo/Rh	○	○	○	●	●	○	●		
Philips	MicroDose	W/Al					●	○		
Planmed	Nuance Classic, Sofie Classic	Mo/Mo	○	●	●	●	●	○	●	
		Mo/Rh	○	○	○	●	●	○	●	
	Nuance, Nuance Excel, Clarity 2D/3D	W/Rh	○	●	○	●	●	○	●	
	W/Ag	○	○		○	○	○	●		
Siemens	Mammomat	Mo/Mo	○	●	●	●	●	○	●	
	Mammomat Inspiration, Mammomat Fusion	Mo/Rh	○	○	○	●	●	○	●	
		W/Rh	○	●	○	●	●	○	●	
	Mammomat Revelation	W/Rh	○	●	○	●	●	○	●	
	W/Ti						○	○		

Note! This list does not include all supported manufacturers or models. Visit <http://www.raysafe.com> to read the product specification.